

For Research Use Only

POLH Polyclonal antibody

Catalog Number: 28133-1-AP



Basic Information

Catalog Number:

28133-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27973

GenBank Accession Number:

BC015742

GeneID (NCBI):

5429

UNIPROT ID:

Q9Y253

Full Name:

polymerase (DNA directed), eta

Calculated MW:

713aa, 78 kDa; 414aa, 46 kDa

Observed MW:

78-85 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : BxPC-3 cells, HeLa cells, COLO 320 cells, HEK-293T cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



BxPC-3 cells were subjected to SDS PAGE followed by western blot with 28133-1-AP (POLH antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.