

For Research Use Only

FOXJ1 Polyclonal antibody

Catalog Number: 28443-1-AP **1 Publications**



Basic Information

Catalog Number: 28443-1-AP	GenBank Accession Number: BC046460	Purification Method: Antigen affinity purification
Size: 900 µg/ml	GeneID (NCBI): 2302	Recommended Dilutions: WB 1:1000-1:5000
Source: Rabbit	UNIPROT ID: Q92949	
Isotype: IgG	Full Name: forkhead box J1	
Immunogen Catalog Number: AG28294	Calculated MW: 45 kDa	
	Observed MW: 45-55 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Caco-2 cells, HeLa cells
Cited Applications: IF	
Species Specificity: Human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xiao-Ling Zou	36929635	Clin Respir J	IF


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28443-1-AP (FOXJ1 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.