

For Research Use Only

RSPO4 Polyclonal antibody

Catalog Number: 28559-1-AP



Basic Information

Catalog Number:

28559-1-AP

Size:

750 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29308

GenBank Accession Number:

NM_001029871

GeneID (NCBI):

343637

UNIPROT ID:

Q2I0M5

Full Name:

R-spondin family, member 4

Calculated MW:

26 kDa

Observed MW:

~19 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : rat brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Rat brain tissue lysates were subjected to SDS PAGE followed by western blot with 28559-1-AP (RSP04 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.