

For Research Use Only

# SOX11 Polyclonal antibody

Catalog Number: 29395-1-AP

Featured Product

1 Publications



## Basic Information

Catalog Number:

29395-1-AP

Size:

1000  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30020

GenBank Accession Number:

NM\_003108

GeneID (NCBI):

6664

UNIPROT ID:

P35716

Full Name:

SRX (sex determining region Y)-box

11

Calculated MW:

47 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : fetal human brain tissue, Y79 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Shaoqin Tu	36910334	Stem Cells Int	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29395-1-AP (SOX11 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.