

For Research Use Only

ZHX3 Polyclonal antibody

Catalog Number: 29397-1-AP

1 Publications



Basic Information

Catalog Number:

29397-1-AP

Size:

1100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30023

GenBank Accession Number:

NM_015035

GeneID (NCBI):

23051

UNIPROT ID:

Q9H4I2

Full Name:

zinc fingers and homeoboxes 3

Calculated MW:

105 kDa

Observed MW:

130 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : HEK-293 cells, Jurkat cells, MCF-7 cells, HeLa cells

Background Information

ZHX3 is a member of the zinc fingers and homeoboxes (ZHX) gene family. ZHX3 is a ubiquitous transcriptional repressor expressed in human cells and tissues.

Notable Publications

Author	Pubmed ID	Journal	Application
Dongwei Yao	39267687	Am J Cancer Res	WB,IF

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29397-1-AP (ZHX3 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.