

For Research Use Only

P54 of ASFV Polyclonal antibody

Catalog Number: 29511-1-AP



Basic Information

Catalog Number: 29511-1-AP	GenBank Accession Number: AAA65354	Purification Method: Antigen affinity purification
Size: 1000 µg/ml	GeneID (NCBI): 22220355	Recommended Dilutions: ELISA 1:556544-1:2226176
Source: Rabbit	UNIPROT ID: Q65194	
Isotype: IgG	Full Name: Envelope protein p54	
Immunogen Catalog Number: AG31444	Calculated MW: 20 kDa	

Applications

Tested Applications: ELISA	Positive Controls: ELISA : FusionProtein,
Species Specificity: Human	

Background Information


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



P54 of ASFV Antibody (29511-1-AP) tested by ELISA. P54 of ASFV(1-65aa) were coated onto microtiter plates at 0.15 µg/well and then incubated with a dilution series of P54 of ASFV Antibody (29511-1-AP). Bound antibodies were detected with HRP conjugated anti-Rabbit IgG followed by incubation with HRP Substrate and then measuring the resulting absorbance at 450 nm.