

For Research Use Only

ICK Polyclonal antibody

Catalog Number: 29728-1-AP



Basic Information

Catalog Number:

29728-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30240

GenBank Accession Number:

BC136420

GeneID (NCBI):

22858

UNIPROT ID:

Q9UPZ9

Full Name:

intestinal cell (MAK-like) kinase

Observed MW:

71 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Species Specificity:

human

Positive Controls:

IP: PC-3 cells,

Background Information

ICK is involved in the control of ciliary length, it may play a key role in the development of multiple organ systems and particularly in cardiac development.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-ICK (IP:29728-1-AP, 4ug;
Detection:29728-1-AP 1:400) with PC-3 cells lysate
1575 ug.