

For Research Use Only

STAP2 Polyclonal antibody

Catalog Number:30516-1-AP



Basic Information

Catalog Number:

30516-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28818

GenBank Accession Number:

BC000795

GeneID (NCBI):

55620

UNIPROT ID:

Q9UGK3

Full Name:

signal transducing adaptor family member 2

Calculated MW:

45 kDa

Observed MW:

50-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HepG2 cells, LNCaP cells, MCF-7 cells, MDA-MB-231 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30516-1-AP (STAP2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.