

For Research Use Only

UVRAG Polyclonal antibody

Catalog Number: 31105-1-AP



Basic Information

Catalog Number:

31105-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34414

GenBank Accession Number:

BC064837

GeneID (NCBI):

7405

UNIPROT ID:

Q9P2Y5

Full Name:

UV radiation resistance associated gene

Calculated MW:

699 aa, 78 kDa

Observed MW:

80-90 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HEK-293 cells, HeLa cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31105-1-AP (UVRAG antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.