

For Research Use Only

RPS11 Polyclonal antibody

Catalog Number: 31204-1-AP



Basic Information

Catalog Number: 31204-1-AP	GenBank Accession Number: BC007945	Purification Method: Antigen affinity Purification
Size: 450 µg/ml	GeneID (NCBI): 6205	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	UNIPROT ID: P62280	
Isotype: IgG	Full Name: ribosomal protein S11	
Immunogen Catalog Number: AG34848	Calculated MW: 158 aa, 18 kDa	
	Observed MW: 18 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells, HeLa cells, MDA-MB-453 cells, mouse testis tissue
Species Specificity: Human, mouse	

Background Information

RPS11 is a ribosomal protein involved in ribosome biogenesis. It is a component of the small ribosomal subunit. The ribosome is a large ribonucleoprotein complex responsible for the synthesis of proteins in the cell.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31204-1-AP (RPS11 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.