

For Research Use Only

LMTK2 Polyclonal antibody

Catalog Number: 31540-1-AP



Basic Information

Catalog Number: 31540-1-AP	GenBank Accession Number: BC131504	Purification Method: Antigen affinity Purification
Size: 400 µg/ml	GeneID (NCBI): 22853	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q8IWU2	
Isotype: IgG	Full Name: Lemur tyrosine kinase 2	
Immunogen Catalog Number: AG35620	Observed MW: 250 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HepG2 cells, U2OS cells, BV-2 cells
Species Specificity: human, mouse	

Background Information

Lemur tyrosine kinase 2 (LMTK2) is a transmembrane protein with serine and threonine but not tyrosine kinase activity.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31540-1-AP (LMTK2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.