

For Research Use Only

TMEM151B Polyclonal antibody

Catalog Number: 31682-1-AP



Basic Information

Catalog Number: 31682-1-AP	GenBank Accession Number: NM_001137560	Purification Method: Antigen affinity Purification
Size: 400 µg/ml	GeneID (NCBI): 441151	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q8IW70	
Isotype: IgG	Full Name: transmembrane protein 151B	
Immunogen Catalog Number: AG34641	Calculated MW: 62 kDa	
	Observed MW: 28 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : K-562 cells, SH-SY5Y cells, TT cells
Species Specificity: human	

Background Information

TMEM151B (Transmembrane protein 151B) may be an integral component of the membrane.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 31682-1-AP (TMEM151B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.