

For Research Use Only

# Jph2 Polyclonal antibody

Catalog Number: 51054-2-AP



## Basic Information

Catalog Number:

51054-2-AP

Size:

433 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0468

GenBank Accession Number:

NM\_021566

GeneID (NCBI):

59091

UNIPROT ID:

Q9ET78

Full Name:

junctophilin 2

Calculated MW:

75 kDa 45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

## Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Positive Controls:

IHC : human pancreas cancer tissue,

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

## Background Information

### Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human pancreas cancer using 51054-2-AP (Jph2 antibody) at dilution of 1:50 (under 10x lens).