

For Research Use Only

PI3K p87 Polyclonal antibody

Catalog Number: 55311-1-AP **1 Publications**



Basic Information

Catalog Number:

55311-1-AP

Size:

510 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001010855

GeneID (NCBI):

146850

UNIPROT ID:

Q5UE93

Full Name:

phosphoinositide-3-kinase,
regulatory subunit 6

Calculated MW:

84 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, mouse

Background Information

PIK3R6, also named as C17orf38 and p87PIKAP, is a regulatory subunit of the PI3K gamma complex. It acts as an adapter to drive activation of PIK3CG by G beta gamma proteins. This antibody is specific to PIK3R5.

Notable Publications

Author	Pubmed ID	Journal	Application
Wenbin Huang	34768109	Neoplasia	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data