

For Research Use Only

TBL1XR1 Polyclonal antibody

Catalog Number: 55312-1-AP **4 Publications**



Basic Information

Catalog Number:

55312-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_024665

GeneID (NCBI):

79718

UNIPROT ID:

Q9BZK7

Full Name:

transducin (beta)-like 1 X-linked receptor 1

Calculated MW:

56 kDa

Observed MW:

50-65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB: mouse brain tissue, MCF-7 cells, rat brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yuan Cheng Li	33102219	Front Oncol	IHC
Ling-Min Liao	30439646	Mol Ther Nucleic Acids	WB, IHC
Yue Zheng	35845072	Front Cardiovasc Med	

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 55312-1-AP (TBL1XR1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.