

For Research Use Only

DHCR24 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60028-2-PBS



Basic Information

Catalog Number: 60028-2-PBS	GenBank Accession Number: BC004375	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 1718	CloneNo.: 3A12F5
Source: Mouse	UNIPROT ID: Q15392	
Isotype: IgG1	Full Name: 24-dehydrocholesterol reductase	
Immunogen Catalog Number: AG0756	Calculated MW: 60 kDa, 516 aa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51020-1, DHCR24 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60028-2-PBS. Detection antibody: 60028-1-PBS. Standard:Ag0756. Range: 0.098-100 ng/mL.