

For Research Use Only

ATF4 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60035-2-PBS



Basic Information

Catalog Number: 60035-2-PBS	GenBank Accession Number: BC022088	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 468	CloneNo.: 3E11C11
Source: Mouse	UNIPROT ID: P18848	
Isotype: IgG2b	Full Name: activating transcription factor 4 (tax-responsive enhancer element B67)	
Immunogen Catalog Number: AG1279	Calculated MW: 39 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51019-2, ATF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60035-2-PBS. Detection antibody: 60035-3-PBS. Standard: Ag16665. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP51019-3, ATF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60035-2-PBS. Detection antibody: 60706-3-PBS. Standard: Ag16665. Range: 0.391-100 ng/mL.