

For Research Use Only

HSD17B10 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60069-4-PBS



Basic Information

Catalog Number:

60069-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG1020

GenBank Accession Number:

BC008708

GeneID (NCBI):

3028

UNIPROT ID:

Q99714

Full Name:

hydroxysteroid (17-beta)
dehydrogenase 10

Calculated MW:

26 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

5F2F6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50654-2, HSD17B10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60069-1-PBS. Detection antibody: 60069-4-PBS. Standard: Ag1020. Range: 0.098-50 ng/mL.