

For Research Use Only

IL-18 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60070-2-PBS



Basic Information

Catalog Number: 60070-2-PBS	GenBank Accession Number: BC007461	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3606	CloneNo.: 5F9B3
Source: Mouse	UNIPROT ID: Q14116	
Isotype: IgG1	Full Name: IL-18	
Immunogen Catalog Number: AG1063	Calculated MW: 193 aa, 22 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP51140-1, IL-18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60070-2-PBS. Detection antibody: 60070-3-PBS. Standard: Ag1063. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP51140-2, IL-18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60070-2-PBS. Detection antibody: 60070-4-PBS. Standard: Ag1063. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP51140-3, IL-18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60070-1-PBS. Detection antibody: 60070-2-PBS. Standard: Ag1063. Range: 0.098-100 ng/mL.