

For Research Use Only

NENF Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60131-2-PBS



Basic Information

Catalog Number: 60131-2-PBS	GenBank Accession Number: BC008823	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 29937	CloneNo.: 5G9F12
Source: Mouse	UNIPROT ID: Q9UMX5	
Isotype: IgG2a	Full Name: neuron derived neurotrophic factor	
Immunogen Catalog Number: AG8387	Calculated MW: 172 aa, 19 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50685-1, NENF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60131-1-PBS. Detection antibody: 60131-2-PBS. Standard:Ag8387. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50685-3, NENF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60131-3-PBS. Detection antibody: 60131-2-PBS. Standard:Ag8387. Range: 0.098-12.5 ng/mL.