

For Research Use Only

MSH2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 60161-3-PBS

Basic Information

Catalog Number: 60161-3-PBS	GenBank Accession Number: BC021566	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4436	CloneNo.: 1F5B10
Source: Mouse	UNIPROT ID: P43246	
Isotype: IgG1	Full Name: mutS homolog 2, colon cancer, nonpolyposis type 1 (E. coli)	
Immunogen Catalog Number: AG7835	Calculated MW: 105 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50235-2, MSH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60161-2-PBS. Detection antibody: 60161-3-PBS. Standard: Ag7835. Range: 0.098-100 ng/mL.