

For Research Use Only

CD122/IL-2RB Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60174-3-PBS



Basic Information

Catalog Number: 60174-3-PBS	GenBank Accession Number: BC025691	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3560	CloneNo.: 3F12D9
Source: Mouse	ENSEMBL Gene ID: ENSG00000100385	
Isotype: IgG1	UNIPROT ID: P14784	
Immunogen Catalog Number: AG3929	Full Name: interleukin 2 receptor, beta	
	Calculated MW: 551 aa, 61 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50451-1, CD122/IL2RB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60174-3-PBS. Detection antibody: 60174-1-PBS. Standard:Ag3929. Range: 0.098-25 ng/mL.



Cytometric bead array standard curve of MP50451-3, CD122/IL2RB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60174-2-PBS. Detection antibody: 60174-3-PBS. Standard:Ag3929. Range: 0.098-3.125 ng/mL.