

For Research Use Only

# CD3D Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60194-3-PBS



## Basic Information

Catalog Number:

60194-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG10207

GenBank Accession Number:

BC070321

GeneID (NCBI):

915

UNIPROT ID:

P04234

Full Name:

CD3d molecule, delta (CD3-TCR complex)

Calculated MW:

171 aa, 19 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

4D9C1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50691-1, CD3D/CD3 Delta Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60194-2-PBS. Detection antibody: 60194-3-PBS. Standard: Ag10207. Range: 0.781-100 ng/mL