

For Research Use Only

Desmin Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60226-2-PBS

Basic Information

Catalog Number:

60226-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9675

GenBank Accession Number:

BC032116

GeneID (NCBI):

1674

UNIPROT ID:

P17661

Full Name:

desmin

Calculated MW:

470 aa, 54 kDa

Purification Method:

Protein G purification

CloneNo.:

6E2G10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50204-1, Desmin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60226-2-PBS. Detection antibody: HRP-conjugated 60226-3-PBS. Standard: ag9675. Range: 0.195-25 ng/mL