

For Research Use Only

# p504S,AMACR Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60240-2-PBS

## Basic Information

Catalog Number:

60240-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8720

GenBank Accession Number:

BC009471

GeneID (NCBI):

23600

UNIPROT ID:

Q9UHK6

Full Name:

alpha-methylacyl-CoA racemase

Calculated MW:

382 aa, 42 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5F4D7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50263-1, p504S, AMACR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60240-2-PBS. Detection antibody: 60240-3-PBS. Standard: Ag8720. Range: 1.563-200 ng/mL.