

For Research Use Only

p504S,AMACR Monoclonal antibody, PBS Only (Detector)



Catalog Number:60240-4-PBS

Basic Information

| | | |
|-------------------------------------|---|---|
| Catalog Number: 60240-4-PBS | GenBank Accession Number: BC009471 | Purification Method: Protein G Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 23600 | CloneNo.: 1C5B12 |
| Source: Mouse | UNIPROT ID: Q9UHK6 | |
| Isotype: IgG1 | Full Name: alpha-methylacyl-CoA racemase | |
| Immunogen Catalog Number: AG8720 | Calculated MW: 382 aa, 42 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50263-2, p504S,AMACR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60240-3-PBS. Detection antibody: 60240-4-PBS. Standard:Ag8720. Range: 1.563-200 ng/mL