

For Research Use Only

OCT4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60242-2-PBS

Basic Information

Catalog Number:

60242-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1794

GenBank Accession Number:

BC020712

GeneID (NCBI):

5460

UNIPROT ID:

Q01860

Full Name:

POU class 5 homeobox 1

Calculated MW:

39 kDa

Purification Method:

Protein A purification

CloneNo.:

4B7F2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50231-1, OCT4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60242-2-PBS. Detection antibody: 60242-3-PBS. Standard:Ag1794. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50231-2, OCT4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60242-2-PBS. Detection antibody: 60242-4-PBS. Standard:Ag1794. Range: 0.098-100 ng/mL.