

For Research Use Only

# OCT4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 60242-4-PBS

## Basic Information

<b>Catalog Number:</b> 60242-4-PBS	<b>GenBank Accession Number:</b> BC020712	<b>Purification Method:</b> Protein A Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 5460	<b>CloneNo.:</b> 1F4C10
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q01860	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> POU class 5 homeobox 1	
<b>Immunogen Catalog Number:</b> AG1794	<b>Calculated MW:</b> 39 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50231-2, OCT4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60242-2-PBS. Detection antibody: 60242-4-PBS. Standard:Ag1794. Range: 0.098-100 ng/mL.