

For Research Use Only

KRT15 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 60247-2-PBS

Basic Information

Catalog Number: 60247-2-PBS	GenBank Accession Number: BC002641	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3866	CloneNo.: 4B7B6
Source: Mouse	UNIPROT ID: P19012	
Isotype: IgG2b	Full Name: keratin 15	
Immunogen Catalog Number: AG0185	Calculated MW: 49 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50002-1, KRT15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60247-1-PBS. Detection antibody: HRP-conjugated 60247-2-PBS. Standard: Ag0185. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP50002-1, KRT15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60247-1-PBS. Detection antibody: 60247-2-PBS. Standard: Ag0185. Range: 1.563-100 ng/mL