

For Research Use Only

HADHA Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60250-4-PBS



Basic Information

Catalog Number:

60250-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG1211

GenBank Accession Number:

BC009235

GeneID (NCBI):

3030

UNIPROT ID:

P40939

Full Name:

hydroxyacyl-Coenzyme A
dehydrogenase/3-ketoacyl-
Coenzyme A thiolase/enoyl-
Coenzyme A hydratase (trifunctional
protein), alpha subunit

Calculated MW:

83 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2D2C10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50412-2, HADHA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60250-3-PBS. Detection antibody: 60250-4-PBS. Standard:Ag1211. Range: 3.125-100 ng/mL.