

For Research Use Only

COXIV Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60251-3-PBS



Basic Information

Catalog Number:

60251-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20551

GenBank Accession Number:

BC021236

GeneID (NCBI):

1327

UNIPROT ID:

P13073

Full Name:

cytochrome c oxidase subunit IV

isoform 1

Calculated MW:

19.6 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5G2G12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50551-1, COXIV Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60251-2-PBS. Detection antibody: 60251-3-PBS. Standard:Ag1640. Range: 0.391-25 ng/mL



Cytometric bead array standard curve of MP50551-2, COXIV Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60251-1-PBS. Detection antibody: 60251-3-PBS. Standard:Ag1640. Range: 0.781-50 ng/mL