

For Research Use Only

COXIV Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60251-4-PBS



Basic Information

Catalog Number: 60251-4-PBS	GenBank Accession Number: BC021236	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1327	CloneNo.: 2F3C3
Source: Mouse	UNIPROT ID: P13073	
Isotype: IgG1	Full Name: cytochrome c oxidase subunit IV isoform 1	
Immunogen Catalog Number: AG20551	Calculated MW: 19.6 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50551-3, COXIV Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66110-1-PBS. Detection antibody: 60251-4-PBS. Standard:Ag1640. Range: 1.563-200 ng/mL.