

For Research Use Only

CD14 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 60253-3-PBS

Basic Information

Catalog Number: 60253-3-PBS	GenBank Accession Number: BC010507	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 929	CloneNo.: 5G12F12
Source: Mouse	ENSEMBL Gene ID: ENSG00000170458	
Isotype: IgG1	UNIPROT ID: P08571	
Immunogen Catalog Number: AG10693	Full Name: CD14 molecule	
	Calculated MW: 375 aa, 40 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50282-1, CD14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60253-2-PBS. Detection antibody: 60253-3-PBS. Standard: Ag10693. Range: 0.781-200 ng/mL.