

For Research Use Only

PTEN Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60300-2-PBS

Basic Information

Catalog Number: 60300-2-PBS	GenBank Accession Number: BC005821	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5728	CloneNo.: 1B1A9
Source: Mouse	UNIPROT ID: P60484	
Isotype: IgG2a	Full Name: phosphatase and tensin homolog	
Immunogen Catalog Number: AG17274	Calculated MW: 47 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50170-1, PTEN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60300-2-PBS. Detection antibody: 60300-3-PBS. Standard: ag17274. Range: 0.391-12.5 ng/mL