

For Research Use Only

c-SRC Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60315-3-PBS



Basic Information

Catalog Number:

60315-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22306

GenBank Accession Number:

BC011566

GeneID (NCBI):

6714

ENSEMBL Gene ID:

ENSG00000197122

UNIPROT ID:

P12931

Full Name:

v-src sarcoma (Schmidt-Ruppin A-2)
viral oncogene homolog (avian)

Calculated MW:

60 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5C1B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51096-1, c-SRC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60315-2-PBS. Detection antibody: 60315-3-PBS. Standard:Ag22306. Range: 3.125-100 ng/mL.