

For Research Use Only

c-SRC Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60315-5-PBS



Basic Information

Catalog Number: 60315-5-PBS	GenBank Accession Number: BC011566	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 6714	CloneNo.: 5C6B3
Source: Mouse	ENSEMBL Gene ID: ENSG00000197122	
Isotype: IgG1	UNIPROT ID: P12931	
Immunogen Catalog Number: AG22306	Full Name: v-src sarcoma (Schmidt-Ruppin A-2) viral oncogene homolog (avian)	
	Calculated MW: 60 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51096-2, c-SRC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60315-4-PBS. Detection antibody: 60315-5-PBS. Standard: Ag22306. Range: 3.125-100 ng/mL.