

For Research Use Only

ALK,CD246 Monoclonal antibody, PBS Only (Capture)



Catalog Number:60321-4-PBS

Basic Information

Catalog Number: 60321-4-PBS	GenBank Accession Number: NM_004304	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 238	CloneNo.: 6G3G6
Source: Mouse	UNIPROT ID: Q9UM73	
Isotype: IgG2a	Full Name: anaplastic lymphoma receptor tyrosine kinase	
Immunogen Catalog Number: AG21493	Calculated MW: 1620 aa, 176 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50317-2, ALK,CD246 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60321-4-PBS. Detection antibody: 60321-3-PBS. Standard:Ag21493. Range: 3.125-100 ng/mL