

For Research Use Only

TMEM106B Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60333-4-PBS

Basic Information

Catalog Number: 60333-4-PBS	GenBank Accession Number: BC033901	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 54664	CloneNo.: 1E12D4
Source: Mouse	UNIPROT ID: Q9NUM4	
Isotype: IgG1	Full Name: transmembrane protein 106B	
Immunogen Catalog Number: AG21448	Calculated MW: 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50326-2, TMEM106B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60333-4-PBS. Detection antibody: 60333-5-PBS. Standard: Ag21448. Range: 0.391-100 ng/mL.