

For Research Use Only

CGB,hCG Monoclonal antibody, PBS Only (Capture)



Catalog Number:60334-2-PBS

Basic Information

Catalog Number: 60334-2-PBS	GenBank Accession Number: BC022796	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1082	CloneNo.: 2F3D4
Source: Mouse	UNIPROT ID: PODN86	
Isotype: IgG2b	Full Name: CG beta	
Immunogen Catalog Number: AG2191	Calculated MW: 165 aa, 18 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50292-1, CGB,hCG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60334-2-PBS. Detection antibody: 60334-1-PBS. Standard:Ag2191. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50292-2, CGB,hCG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60334-2-PBS. Detection antibody: 60334-3-PBS. Standard:Ag2191. Range: 0.098-100 ng/mL.