

For Research Use Only

# E-cadherin Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60335-2-PBS

## Basic Information

Catalog Number:

60335-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG15085

GenBank Accession Number:

BC141838

GeneID (NCBI):

999

UNIPROT ID:

P12830

Full Name:

cadherin 1, type 1, E-cadherin (epithelial)

Calculated MW:

882 aa, 97 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2H5A10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50315-1, E-cadherin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60335-2-PBS. Detection antibody: 60335-3-PBS. Standard: Ag15085. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP50315-2, E-cadherin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60335-2-PBS. Detection antibody: 60335-4-PBS. Standard: Ag15085. Range: 0.195-12.5 ng/mL