

For Research Use Only

E-cadherin Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 60335-4-PBS

Basic Information

Catalog Number:

60335-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15085

GenBank Accession Number:

BC141838

GeneID (NCBI):

999

UNIPROT ID:

P12830

Full Name:

cadherin 1, type 1, E-cadherin (epithelial)

Calculated MW:

882 aa, 97 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E4E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50315-2, E-cadherin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60335-2-PBS. Detection antibody: 60335-4-PBS. Standard: Ag15085. Range: 0.195-12.5 ng/mL



Cytometric bead array standard curve of MP50315-3, E-cadherin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60335-4-PBS. Detection antibody: 60335-5-PBS. Standard: Ag15085. Range: 1.563-50 ng/mL