

For Research Use Only

# BSAP, PAX5 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60349-4-PBS

## Basic Information

Catalog Number:

60349-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16023

GenBank Accession Number:

BC156927

GeneID (NCBI):

5079

UNIPROT ID:

Q02548

Full Name:

paired box 5

Calculated MW:

391 aa, 42 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3D2H5

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50228-3, BSAP,PAX5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60349-4-PBS. Detection antibody: 60349-2-PBS. Standard:Ag16023. Range: 0.781-100 ng/mL