

For Research Use Only

MYO5B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60402-1-PBS



Basic Information

Catalog Number: 60402-1-PBS	GenBank Accession Number: NM_001080467	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4645	CloneNo.: 1D1E2
Source: Mouse	UNIPROT ID: Q9ULV0	
Isotype: IgG2a	Full Name: myosin VB	
Immunogen Catalog Number: AG32670	Calculated MW: 214KD	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50502-1, MYO5B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60402-1-PBS. Detection antibody: 60402-2-PBS. Standard: Ag32670. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50502-2, MYO5B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60402-1-PBS. Detection antibody: 60402-3-PBS. Standard: Ag32670. Range: 1.563-100 ng/mL.