

For Research Use Only

# MYO5B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60402-7-PBS



## Basic Information

|   |  |  |
|---|--|--|
| <b>Catalog Number:</b><br>60402-7-PBS       | <b>GenBank Accession Number:</b><br>NM_001080467 | <b>Purification Method:</b><br>Protein A Magarose purification |
| <b>Size:</b><br>1 mg/ml                     | <b>GeneID (NCBI):</b><br>4645                    | <b>CloneNo.:</b><br>2A9F9                                      |
| <b>Source:</b><br>Mouse                     | <b>UNIPROT ID:</b><br>Q9ULV0                     |  |
| <b>Isotype:</b><br>IgG2a                    | <b>Full Name:</b><br>myosin VB                   |  |
| <b>Immunogen Catalog Number:</b><br>AG32670 | <b>Calculated MW:</b><br>214KD                   |  |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50502-4, MYO5B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60402-6-PBS. Detection antibody: 60402-7-PBS. Standard:Ag32670. Range: 1.563-100 ng/mL.