

For Research Use Only

FGFR2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60403-1-PBS



Basic Information

Catalog Number:

60403-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

2263

UNIPROT ID:

P21802-3

Full Name:

fibroblast growth factor receptor 2

Calculated MW:

92kd

Purification Method:

Protein A Magarose purification

CloneNo.:

3G11A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50495-1, FGFR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60403-1-PBS. Detection antibody: 60403-2-PBS. Standard:Eg0175. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50495-3, FGFR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60403-1-PBS. Detection antibody: 60403-5-PBS. Standard:Eg0175. Range: 0.391-25 ng/mL.