

For Research Use Only

USP19 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60404-1-PBS



Basic Information

Catalog Number:

60404-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG22645

GenBank Accession Number:

BC146752

GeneID (NCBI):

10869

UNIPROT ID:

O94966

Full Name:

ubiquitin specific peptidase 19

Calculated MW:

1318 aa, 146 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F2F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50506-1, USP19 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60404-1-PBS. Detection antibody: 60404-2-PBS. Standard: Ag22645. Range: 0.391-6.25 ng/mL