

For Research Use Only

ALDH1L2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60405-2-PBS



Basic Information

Catalog Number:

60405-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG17340

GenBank Accession Number:

BC103935

GeneID (NCBI):

160428

UNIPROT ID:

Q3SY69

Full Name:

aldehyde dehydrogenase 1 family, member L2

Calculated MW:

923 aa, 102 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F3H1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50511-1, ALDH1L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60405-1-PBS. Detection antibody: 60405-2-PBS. Standard:Ag17340. Range: 0.391-12.5 ng/mL



Cytometric bead array standard curve of MP50511-2, ALDH1L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60405-3-PBS. Detection antibody: 60405-2-PBS. Standard:Ag17340. Range: 0.391-25 ng/mL