

For Research Use Only

CHAC1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60406-1-PBS



Basic Information

Catalog Number:

60406-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7527

GenBank Accession Number:

BC001847

GeneID (NCBI):

79094

UNIPROT ID:

Q9BUX1

Full Name:

ChaC, cation transport regulator homolog 1 (E. coli)

Calculated MW:

24 kDa

Purification Method:

Protein G purification

CloneNo.:

2H1D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50512-1, CHAC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60406-1-PBS. Detection antibody: 60406-2-PBS. Standard:Ag7527. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50512-2, CHAC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60406-1-PBS. Detection antibody: 60406-3-PBS. Standard:Ag7527. Range: 0.195-50 ng/mL.