

For Research Use Only

SLC44A1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60408-1-PBS



Basic Information

Catalog Number: 60408-1-PBS	GenBank Accession Number: BC049203	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 23446	CloneNo.: 1B2C4
Source: Mouse	UNIPROT ID: Q8WWI5	
Isotype: IgG2a	Full Name: solute carrier family 44, member 1	
Immunogen Catalog Number: AG6024	Calculated MW: 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50519-1, SLC44A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60408-1-PBS. Detection antibody: 60408-2-PBS. Standard: Ag6024. Range: 0.391-100 ng/mL.