

For Research Use Only

PIP4K2A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60415-1-PBS



Basic Information

Catalog Number:

60415-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG34702

GenBank Accession Number:

BC018034

GeneID (NCBI):

5305

UNIPROT ID:

P48426

Full Name:

phosphatidylinositol-5-phosphate 4-kinase, type II, alpha

Calculated MW:

406 aa, 46 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3H2A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50530-1, PIP4K2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60415-1-PBS. Detection antibody: 60415-2-PBS. Standard: Ag34702. Range: 0.391-50 ng/mL